

## Product Datasheet

### Rat BCL2L1 ELISA Kit (orb410607)

|                          |                                                   |
|--------------------------|---------------------------------------------------|
| <b>Description</b>       | Rat Bcl-2-like protein 1(Bcl2-L-1/Bcl-X)ELISA Kit |
| <b>Reactivity</b>        | Rat                                               |
| <b>Range</b>             | 46.88 pg/mL-3000 pg/mL                            |
| <b>Target</b>            | BCL2L1                                            |
| <b>Note</b>              | For research use only                             |
| <b>Application notes</b> | Detection Wavelength: 450 nm                      |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates, cell lysates   |
| <b>Uniprot ID</b>        | <b>P53563</b>                                     |
| <b>Sensitivity</b>       | 11.72 pg/mL                                       |
| <b>Expiration Date</b>   | Please enquire.                                   |