

Product Datasheet

Rat CCL6 ELISA Kit (orb410615)

Description	Rat C-C motif chemokine 6 (CCL6) ELISA kit
Reactivity	Rat
Range	0.156 ng/mL-10 ng/mL
Target	CCL6
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q68FP3
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.