

Product Datasheet

Rat Cluster of Differentiation 68 ELISA Kit (orb410622)

Description	Rat Cluster of Differentiation 68, CD68 ELISA kit
Reactivity	Rat
Range	0.625 ng/mL-40 ng/mL
Target	CD68
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q4FZY1
Sensitivity	0.156 ng/mL
Expiration Date	Please enquire.