

Product Datasheet

Rat CXCR4 ELISA Kit (orb410625)

Description	Rat CXC-chemokine receptor 4(CXCR4)ELISA Kit
Reactivity	Rat
Range	0.156 ng/mL-10 ng/mL
Target	CXCR4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	O08565
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.