

Product Datasheet

Rat CRP ELISA Kit (orb410626)

Description	Rat C-Reactive Protein,CRP ELISA Kit
Reactivity	Rat
Range	31.25 ng/mL-2000 ng/mL
Target	CRP
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, cell lysates
Uniprot ID	P48199
Sensitivity	7.81 ng/mL
Expiration Date	Please enquire.