

Product Datasheet

Rat Flt-3L ELISA Kit (orb410635)

Description	Rat FMS-like tyrosine kinase 3 ligand, Flt-3L ELISA Kit
Reactivity	Rat
Range	15.6 pg/mL-1000 pg/mL
Target	Flt-3L
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q91ZT1
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.