

Product Datasheet

Rat NAG ELISA Kit (orb410651)

Description	Rat N-acetyl- β -D-glucosaminidase,NAG ELISA Kit
Reactivity	Rat
Range	1.56 mIU/mL-100 mIU/mL
Target	NAG
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, urine
Uniprot ID	Q6AXR4
Sensitivity	0.39 mIU/mL
Expiration Date	Please enquire.