

Product Datasheet

Rat RUNX2 ELISA Kit (orb410659)

Description	Rat Runt-related transcription factor 2(RUNX2) ELISA kit
Reactivity	Rat
Range	3.12 ng/mL-200 ng/mL
Target	RUNX2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9Z2J9
Sensitivity	0.78 ng/mL
Expiration Date	Please enquire.