

Product Datasheet

Rat ORM1 ELISA Kit (orb410675)

Description	Rat α 1-Acid glycoprotein(α 1-AGP) ELISA Kit
Reactivity	Rat
Range	3.125 ng/mL-800 ng/mL
Target	ORM1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, urine, cell culture supernates, tissue homogenates
Uniprot ID	P02764
Sensitivity	3.125 ng/mL
Expiration Date	Please enquire.