

## Product Datasheet

### Rat IFN $\alpha$ ELISA Kit (orb410679)

|                          |                                                               |
|--------------------------|---------------------------------------------------------------|
| <b>Description</b>       | Rat Interferon $\alpha$ , IFN- $\alpha$ ELISA Kit             |
| <b>Reactivity</b>        | Rat                                                           |
| <b>Range</b>             | 31.25 pg/mL-2000 pg/mL                                        |
| <b>Target</b>            | IFNA1                                                         |
| <b>Note</b>              | For research use only                                         |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                  |
| <b>Sample Types</b>      | serum, plasma, cell culture supernates,<br>tissue homogenates |
| <b>Uniprot ID</b>        | <b>P05011</b>                                                 |
| <b>Sensitivity</b>       | 7.81 pg/mL                                                    |
| <b>Expiration Date</b>   | Please enquire.                                               |