

Product Datasheet

Rat amyloid beta peptide 1-42, A β 1-42 ELISA Kit (orb410690)

Description	Rat amyloid beta peptide 1-42, A β 1-42 ELISA Kit
Reactivity	Rat
Range	0.9 pg/mL-60 pg/mL
Target	A β 1-42
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Sensitivity	0.225 pg/mL
Expiration Date	Please enquire.