

Product Datasheet

Rat GGT1 ELISA Kit (orb410699)

Description	Rat Gamma-glutamyltranspeptidase 1 (GGT1) ELISA kit
Reactivity	Rat
Range	4.7 mIU/mL-300 mIU/mL
Target	GGT1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P07314
Sensitivity	1.2 mIU/mL
Expiration Date	Please enquire.