

Product Datasheet

Rat ST2 ELISA Kit (orb410721)

Description	Rat IL-1 Receptor Like 1(IL1RL1)ELISA Kit
Reactivity	Rat
Range	15.6 pg/mL-1000 pg/mL
Target	IL1RL1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q62611
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.