

Product Datasheet

Rat SOD3 ELISA Kit (orb410749)

Description	Rat extracellular superoxide dismutase [Cu-Zn] (Cu/Zn-SOD/SOD3)ELISA kit
Reactivity	Rat
Range	15.6 pg/mL-1000 pg/mL
Target	SOD3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q08420
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.