

Product Datasheet

Rat SOD2 ELISA Kit (orb410768)

Description	Rat Superoxide dismutase [Mn], mitochondrial(SOD2) ELISA kit
Reactivity	Rat
Range	1.87 pg/mL-120 pg/mL
Target	SOD2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P07895
Sensitivity	0.47 pg/mL
Expiration Date	Please enquire.