

Product Datasheet

Rat FGF23 ELISA Kit (orb410773)

Description	Rat fibroblast growth factor 23 (FGF23) ELISA Kit
Reactivity	Rat
Range	1.56 pg/mL-100 pg/mL
Target	FGF23
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q8VI82
Sensitivity	0.39 pg/mL
Expiration Date	Please enquire.