

Product Datasheet

Rat rudimental bovine serum albumin(BSA) check-up ELISA kit [Out of stock] (orb410805)

Description	Rat rudimental bovine serum albumin check-up ELISA Kit
Reactivity	Rat
Range	0.312 µg/mL-80 µg/mL
Target	rudimental bovine serum albumin check-up
Note	For research use only
Application notes	• Sample Type: serum, plasma,tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Sensitivity	0.312 µg/mL
Expiration Date	Please enquire.