

Product Datasheet

Rat SFTPA1 ELISA Kit (orb410833)

Description	Rat Pulmonary surfactant-associated protein A,SP-A ELISA Kit
Reactivity	Rat
Range	31.25 pg/ml - 2000 pg/ml
Target	SFTPA1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P08427
Sensitivity	7.8 pg/ml
Expiration Date	Please enquire.