

Product Datasheet

Rat GAD autoantibody (IgG) ELISA Kit (orb410866)

Description	Rat glutamic acid decarboxylase(GAD) autoantibody(IgG) ELISA Kit
Reactivity	Rat
Range	7.8 pg/mL-500 pg/mL
Target	GAD auto-Ab (IgG)
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma
Sensitivity	1.95 pg/mL
Expiration Date	Please enquire.