

Product Datasheet

Rat SREBP-1C ELISA Kit (orb410887)

Description	Rat Sterol regulatory element-binding protein 1C(SREBP-1C) ELISA Kit
Reactivity	Rat
Range	62.5 pg/mL-4000 pg/mL
Target	SREBP-1C
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P56720
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.