

Product Datasheet

SFI1 antibody (orb41089)

Description	SFI1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Protein SFI1 homolog protein (1003-1242AA)
Target	SFI1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A8K8P3
Dilution Range	IHC-P:1:20-200
Expiration Date	12 months from date of receipt.