

## Product Datasheet

### Rat FRA ELISA Kit (orb410897)

|                          |                                   |
|--------------------------|-----------------------------------|
| <b>Description</b>       | Rat fructosamin (FRA) ELISA kit   |
| <b>Reactivity</b>        | Rat                               |
| <b>Range</b>             | 31.25 nmol/mL-2000 nmol/mL        |
| <b>Target</b>            | FRA                               |
| <b>Note</b>              | For research use only             |
| <b>Application notes</b> | Detection Wavelength: 450 nm      |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates |
| <b>Sensitivity</b>       | 7.81 nmol/mL                      |
| <b>Expiration Date</b>   | Please enquire.                   |