

Product Datasheet

Rat PPARG ELISA Kit (orb410902)

Description	Rat Peroxisome Proliferator-activated receptor gamma;;PPAR-gamma; ELISA Kit
Reactivity	Rat
Range	0.156 ng/mL-10 ng/mL
Target	PPARG
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O88275
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.