

## Product Datasheet

### Rat EPO ELISA Kit (orb410908)

|                          |                                      |
|--------------------------|--------------------------------------|
| <b>Description</b>       | Rat Erythropoietin,EPO ELISA Kit     |
| <b>Reactivity</b>        | Rat                                  |
| <b>Range</b>             | 0.156 ng/mL-10 ng/mL                 |
| <b>Target</b>            | EPO                                  |
| <b>Note</b>              | For research use only                |
| <b>Application notes</b> | Detection Wavelength: 450 nm         |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates |
| <b>Uniprot ID</b>        | <b>P29676</b>                        |
| <b>Sensitivity</b>       | 0.039 ng/mL                          |
| <b>Expiration Date</b>   | Please enquire.                      |