

## Product Datasheet

### Rat GHRL ELISA Kit (orb410922)

|                          |                                        |
|--------------------------|----------------------------------------|
| <b>Description</b>       | Rat ghrelin (GHRL) ELISA kit           |
| <b>Reactivity</b>        | Rat                                    |
| <b>Range</b>             | 7.8 pg/mL-500 pg/mL                    |
| <b>Target</b>            | GHRL                                   |
| <b>Note</b>              | For research use only                  |
| <b>Application notes</b> | Detection Wavelength: 450 nm           |
| <b>Sample Types</b>      | serum, plasma, cell culture supernates |
| <b>Uniprot ID</b>        | <b>Q9QYH7</b>                          |
| <b>Sensitivity</b>       | 1.95 pg/mL                             |
| <b>Expiration Date</b>   | Please enquire.                        |