

Product Datasheet

Rat MMP1 ELISA Kit (orb410930)

Description	Rat matrix metalloproteinase 1(MMP-1/Collagenase I) ELISA kit
Reactivity	Rat
Range	0.235 ng/ml - 15 ng/ml
Target	MMP1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P81563
Sensitivity	0.058 ng/ml
Expiration Date	Please enquire.