

Product Datasheet

Rat MMP2 ELISA Kit (orb410934)

Description	Rat matrix metalloproteinase 2/Gelatinase A,MMP-2 ELISA kit
Reactivity	Rat
Range	0.625 ng/mL-40 ng/mL
Target	MMP2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P33436
Sensitivity	0.156 ng/mL
Expiration Date	Please enquire.