

Product Datasheet

Rat MMP8 ELISA Kit (orb410938)

Description	Rat matrix metallopeptidase 8 (neutrophil collagenase) (MMP8) ELISA kit
Reactivity	Rat
Range	0.312 ng/mL-20 ng/mL
Target	MMP8
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	O88766
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.