

Product Datasheet

Rat MMP9 ELISA Kit (orb410939)

Description	Rat matrix metalloproteinase 9/Gelatinase B,MMP-9 ELISA Kit
Reactivity	Rat
Range	0.78 ng/ml-50 ng/ml.
Target	MMP9
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates.
Uniprot ID	P50282
Sensitivity	0.195 ng/ml.
Expiration Date	Please enquire.