

## Product Datasheet

### Rat SDF1a/CXCL12 ELISA Kit (orb410942)

|                          |                                                            |
|--------------------------|------------------------------------------------------------|
| <b>Description</b>       | Rat Stromal cell derived factor 1a,SDF-1a ELISA Kit        |
| <b>Reactivity</b>        | Rat                                                        |
| <b>Range</b>             | 0.78 pg/mL-50 pg/mL                                        |
| <b>Target</b>            | SDF1a/CXCL12                                               |
| <b>Note</b>              | For research use only                                      |
| <b>Application notes</b> | Detection Wavelength: 450 nm                               |
| <b>Sample Types</b>      | serum, plasma, cell culture supernates, tissue homogenates |
| <b>Uniprot ID</b>        | <b>Q9QZD1</b>                                              |
| <b>Sensitivity</b>       | 0.195 pg/mL                                                |
| <b>Expiration Date</b>   | Please enquire.                                            |