

Product Datasheet

Rat Haptoglobin ELISA Kit (orb410963)

Description	Rat Haptoglobin,Hpt/HP ELISA KIT
Reactivity	Rat
Range	0.156 µg/mL-10 µg/mL
Target	HP
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates,bronchoalveolar lavage fluid, urine
Uniprot ID	P06866
Sensitivity	0.039 µg/mL
Expiration Date	Please enquire.