



Product Datasheet

UQCC2 antibody (orb41097)


Description	UQCC2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Ubiquinol-cytochrome-c reductase complex assembly factor 2 protein (14-126AA)
Target	UQCC2
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BRT2
Dilution Range	WB:1:1000-5000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human heart tissue.



Immunohistochemical staining of human liver tissue.



Western blot analysis of A431 whole cell lysate.