

Product Datasheet

Rat CSF1 ELISA Kit (orb410973)

Description	Rat Macrophage Colony-Stimulating Factor,M-CSF ELISA kit
Reactivity	Rat
Range	15.6 pg/mL-1000 pg/mL
Target	CSF1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q8JZQ0
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.