

Product Datasheet

Rat CD40LG ELISA Kit (orb410997)

Description	Rat Soluble Cluster of differentiation 40 ligand,sCD40L ELISA Kit
Reactivity	Rat
Range	1.87 ng/mL-120 ng/mL
Target	CD40LG
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9Z2V2
Sensitivity	0.47 ng/mL
Expiration Date	Please enquire.