

Product Datasheet

Rat B72/sCD86 ELISA Kit (orb410999)

Description	Rat soluble Cluster of differentiation 86,sCD86 ELISA Kit
Reactivity	Rat
Range	78 pg/mL-5000 pg/mL
Target	B72/sCD86
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma,tissue homogenates
Uniprot ID	O35531
Sensitivity	19.5 pg/mL
Expiration Date	Please enquire.