

Product Datasheet

RPP25L antibody (orb41100)

Description	RPP25L antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Ribonuclease P protein subunit p25-like protein (1-163AA)
Target	RPP25L
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N5L8
Dilution Range	IHC-P:1:20-200
Expiration Date	12 months from date of receipt.