

Product Datasheet

Rat RANKL ELISA Kit (orb411003)

Description	Rat soluble receptor activator of nuclear factor-kB ligand,sRANKL ELISA Kit
Reactivity	Rat
Range	62.5 pg/mL-4000 pg/mL
Target	TNFSF11
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates
Uniprot ID	Q9ESE2
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.