

Product Datasheet

Rat SRD5A1 ELISA Kit (orb411014)

Description	Rat Steroid 5-Alpha-Reductase 1 (SRD5A1) ELISA Kit
Reactivity	Rat
Range	15.6 pg/mL-1000 pg/mL
Target	SRD5A1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P24008
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.