

## Product Datasheet

### Rat G-CSF ELISA Kit (orb411015)

|                          |                                                            |
|--------------------------|------------------------------------------------------------|
| <b>Description</b>       | Rat Granulocyte Colony Stimulating Factor,G-CSF ELISA Kit  |
| <b>Reactivity</b>        | Rat                                                        |
| <b>Range</b>             | 31.25 pg/mL-2000 pg/mL                                     |
| <b>Target</b>            | G-CSF                                                      |
| <b>Note</b>              | For research use only                                      |
| <b>Application notes</b> | Detection Wavelength: 450 nm                               |
| <b>Sample Types</b>      | serum, plasma, cell culture supernates, tissue homogenates |
| <b>Uniprot ID</b>        | <b>P48750</b>                                              |
| <b>Sensitivity</b>       | 7.81 pg/mL                                                 |
| <b>Expiration Date</b>   | Please enquire.                                            |