

## Product Datasheet

### Rat CSF2 ELISA Kit (orb411016)

|                          |                                                                       |
|--------------------------|-----------------------------------------------------------------------|
| <b>Description</b>       | Rat Granulocyte-Macrophage Colony Stimulating Factor,GM-CSF ELISA Kit |
| <b>Reactivity</b>        | Rat                                                                   |
| <b>Range</b>             | 15.6 pg/ml - 1000 pg/ml                                               |
| <b>Target</b>            | CSF2                                                                  |
| <b>Note</b>              | For research use only                                                 |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                          |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                                     |
| <b>Uniprot ID</b>        | <b>P48750</b>                                                         |
| <b>Sensitivity</b>       | 3.9 pg/ml                                                             |
| <b>Expiration Date</b>   | Please enquire.                                                       |