

Product Datasheet

Rat HSPG2 ELISA Kit (orb411024)

Description	Rat heparan sulfate proteoglycan 2 (HSPG2) ELISA kit
Reactivity	Rat
Range	15.6 pg/mL-1000 pg/mL
Target	HSPG2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, bronchoalveolar lavage fluid, urine, tissue homogenates
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.