

Product Datasheet

Rat SERPINA12 ELISA Kit (orb411057)

Description	Rat Visceral adipose-specific serine protease inhibitor,vaspin ELISA Kit
Reactivity	Rat
Range	31.25 pg/mL-2000 pg/mL
Target	SERPINA12
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q8R4Z1
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.