

Product Datasheet

Rat UGT1A1 ELISA Kit (orb411062)

Description	Rat UDP-glucuronosyltransferase 1(UGT1) ELISA Kit
Reactivity	Rat
Range	31.2 pg/mL-2000 pg/mL
Target	UGT1A1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q64550
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.