

Product Datasheet

Rat PLAUR ELISA Kit (orb411064)

Description	Rat Urokinase-type Plasminogen Activator Receptor, PLAUR/uPAR ELISA kit
Reactivity	Rat
Range	0.156 ng/mL-10 ng/mL
Target	PLAUR
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P49616
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.