

Product Datasheet

Rat SLC2A4 ELISA Kit (orb411083)

Description	Rat glucose transporter 4(GLUT4) ELISA Kit
Reactivity	Rat
Range	0.156 ng/mL-10 ng/mL
Target	SLC2A4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P19357
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.