

Product Datasheet

Rat PCSK9 ELISA Kit (orb411085)

Description	Rat Proprotein convertase subtilisin/kexin type 9(PCSK9) ELISA kit
Reactivity	Rat
Range	93.75 pg/mL-6000 pg/mL
Target	PCSK9
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P59996
Sensitivity	46.88 pg/mL
Expiration Date	Please enquire.