

Product Datasheet

Rat HSP-90 ELISA Kit (orb411107)

Description	Rat Heat shock protein HSP 90-alpha (Hsp90aa1/Hsp86/Hspca) ELISA kit
Reactivity	Rat
Range	0.312 ng/mL-20 ng/mL
Target	Hsp90aa1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P82995
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.