

Product Datasheet

Rat TPH1 ELISA Kit (orb411118)

Description	Rat Tryptophan Hydroxylase (TPH) ELISA Kit
Reactivity	Rat
Range	37.5 pg/mL-2400 pg/mL
Target	TPH1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P09810
Sensitivity	9.4 pg/mL
Expiration Date	Please enquire.