

Product Datasheet

Rat CXCL5 ELISA Kit (orb411121)

Description	Rat Epithelial neutrophil activating peptide 78,ENA-78 ELISA Kit
Reactivity	Rat
Range	0.312 ng/mL-20 ng/mL
Target	CXCL5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P97885
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.