

Product Datasheet

Rat NOS1 ELISA Kit (orb411136)

Description	Rat Neuronal Nitric Oxide Synthase (nNOS/NOS1) ELISA Kit
Reactivity	Rat
Range	31.25 mIU/mL-2000 mIU/mL
Target	NOS1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, cell lysates
Uniprot ID	P29476
Sensitivity	7.81 mLU/mL
Expiration Date	Please enquire.